Integra®

Matriz de Regeneração Dérmica

Sistema de Regeneração
Dérmica

INTEGRA®
LIMIT UNCERTAINTY

Resultados

Caso 1

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Silver FH, Yannas IV, Salzman EW. Glycosaminoglycan inhibition of collagen induced platelet aggregation. *Thromb Res.* 1978 Aug;13(2):267-77.

Silver FH, Yannas IV, Salzman EW. In vitro blood compatibility of glycosaminoglycan-precipitated collagens. *J Biomed Mater Res.* 1979 Sep;13(5):701-16.

Yannas IV, Burke JF. Design of an artificial skin. I. Basic design principles. *J Biomed Mater Res.* 1980 Jan;14(1):65-81.

Yannas IV, Burke JF, Gordon PL, Huang C, Rubenstein RH. Design of an artificial skin. II. Control of chemical composition. *J Biomed Mater Res.* 1980 Mar;14(2):107-32.

Dagalakis N, Flink J, Stasikelis P, Burke JF, Yannas IV. Design of an artificial skin. Part III. Control of pore structure. *J Biomed Mater Res.* 1980 Jul;14(4):511-28.

Yannas IV, Burke JF, Orgill DP, Skrabut EM. Wound tissue can utilize a polymeric template to synthesize a functional extension of skin. *Science.* 1982;215(4529):174-6.

Yannas IV, Lee E, Orgill DP, Skrabut EM, Murphy GF. Synthesis and characterization of a model extracellular matrix that induces partial regeneration of adult mammalian skin. *Proc Natl Acad Sci U.S.A.* 1989 Feb;86(3):933-7.

Orgill DP, Butler CE, Regan JF. Behavior of collagen-GAG matrices as dermal replacement in rodent and porcine models. *Wounds.* 1996;8(5):151-57.

King WW, Lam PK, Liew CT, Ho WS, Li AK. Evaluation of artificial skin (Integra) in a rodent model. *Burns.* 1997;23 Suppl 1(-HD-):S30-2.

Yannas IV. Studies on the biological activity of the dermal regeneration template. *Wound Repair Regeneration.* 1998;6:518-524.

Compton CC, Butler CE, Yannas IV, Warland G, Orgill DP. Organized skin structure is regenerated in vivo from collagen-GAG matrices seeded with autologous keratinocytes. *J Invest Dermatol.* 1998 Jun;110(6):908-16.

Butler CE, Orgill DP, Yannas IV, Compton CC. Effect of keratinocyte seeding of collagen-glycosaminoglycan membranes on the regeneration of skin in a porcine model. *Plast Reconstr Surg.* 1998 May;101(6):1572-9.

Orgill DP, Butler C, Regan JF, Barlow MS, Yannas IV, Compton CC. Vascularized collagen-glycosaminoglycan matrix provides a dermal substrate and improves take of cultured epithelial autografts. *Plast Reconstr Surg.* 1998 Aug;102(2):423-9.

Orgill DP, Yannas IV. Design of an artificial skin. IV. Use of island graft to isolate organ regeneration from scar synthesis and other processes leading to skin wound closure. *J Biomed Mater Res.* 1998;39(4):531-5.

Butler CE, Yannas IV, Compton CC, Correia CA, Orgill DP. Comparison of cultured and uncultured keratinocytes seeded into a collagen-GAG matrix for skin replacements. *Br J Plast Surg.* 1999 Mar;52(2):127-32.

Kremer M, Lang E, Berger AC. Evaluation of dermal-epidermal skin equivalents ("composite-skin") of human keratinocytes in a collagen-glycosaminoglycan matrix (Integra artificial skin). *Br J Plast Surg.* 2000 Sep;53(6):459-65.

Shermak MA, Wong L, Inoue N, Nicol T. Reconstruction of complex cranial wounds with demineralized bone matrix and bilayer artificial skin. *J Craniofac Surg.* 2000 May;11(3):224-31.

Ojeh NO, Frame JD, Navsaria HA. In Vitro Characterization of an Artificial Dermal Scaffold. *Tissue Eng.* 2001 Aug;7(4):457-72.

Grant I, Green C, Martin R. Strategies to improve the take of commercially available collagen/glycosaminoglycan wound repair material investigated in an animal model. *Burns.* 2001 Nov;27(7):699-707.

Lam PK, Chan ES, Liew CT, Lau C, Yen SC, King WW. Combination of a new composite biocompatible skin graft on the neodermis of artificial skin in an animal model. *ANZ J Surg.* 2002 May;72(5):360-3.

Chu CS, McManus AT, Matyevich NP, Goodwin CW, Pruitt BA Jr. Integra as a dermal replacement in a meshed composite skin graft in a rat model: a one-step operative procedure. *J Trauma.* 2002 Jan;52(1):122-9.

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Burke JF, Yannas IV, Quinby WC Jr, Bondoc CC, Jung WK. Successful use of a physiologically acceptable artificial skin in the treatment of extensive burn injury. *Ann Surg.* 1981;194(4):413-28.

Burke JF. Observations on the development and clinical use of artificial skin—an attempt to employ regeneration rather than scar formation in wound healing. *Jpn J Surg.* 1987;17(6):431-8.

Heimbach D, Luterman A, Burke J, Cram A, Herndon D, Hunt J, Jordan M, McManus W, Solem L, Warden G, Zawacki B. Artificial dermis for major burns. A multi-center randomized clinical trial. *Ann Surg.* 1988;208(3):313-20.

Tompkins RG, Hilton JF, Burke JF, Schoenfeld DA, Hegarty MT, Bondoc CC, Quinby WC Jr, Behringer GE, Ackroyd FW. Increased survival after massive thermal injuries in adults: preliminary report using artificial skin. *Crit Care Med.* 1989;17(8):734-40.

Stern R, McPherson M, Longaker MT. Histologic study of artificial skin used in the treatment of full-thickness thermal injury. *J Burn Care Rehabil.* 1990 Jan-Feb;11:1,7-13.

Michaeli D, McPherson M. Immunologic study of thermal injuries. *J Burn Care Rehabil.* 1990 Jan-Feb;11:1,21-6.

Sheridan RL, Hegarty M, Tompkins RG, Burke JF. Artificial skin in massive burns—results to ten years. *Eur J Plast Surg.* 1994;17:91-93.

Lorenz C, Petravic A, Hohl HP, Wessel L, Waag KL. Early wound closure and early reconstruction. Experience with a dermal substitute in a child with 60 percent surface area burn. *Burns.* 1997 Sep;23: 6,505-8.

Besner GE, Klamar JE. Integra Artificial Skin as a useful adjunct in the treatment of purpura fulminans. *J Burn Care Rehabil.* 1998 Jul-Aug;19:4,324-9.

Clayton MC, Bishop JF. Perioperative and postoperative dressing techniques for Integra Artificial Skin: views from two medical centers. *J Burn Care Rehabil.* 1998;19:4,358-63.

Pandya AN, Woodward B, Parkhouse N. The use of cultured autologous keratinocytes with Integra in the resurfacing of acute burns. *Plast Reconstr Surg.* 1998 Sep;102:3,825-8; discussion 829-30.

Ceddi C, Hartmann B, Schepler H, Raff T, Germann G. Grafting of deeply burned problem zones in the lower extremity with a dermal substitute. *Eur J Plast Surg.* 1999;22:119-124.

Boyce ST, Kagan RJ, Meyer NA, Yakuboff KP, Warden GD. The 1999 clinical research award. Cultured skin substitutes combined with Integra Artificial Skin to replace native skin autograft and allograft for the closure of excised full-thickness burns. *J Burn Care Rehabil.* 1999 Nov-Dec;20:6,453-61.

Hunt JA, Moisisdis E, Haertsch P. Initial experience of Integra in the treatment of post-burn anterior cervical neck contracture. *British Journal Plastic Surgery.* 2000 Dec;53(8):652-658.

Berger A, Tanzella U, Machens HG, Liebau J. Administration of Integra on primary burn wounds and unstable secondary scars. *Chirurg.* 2000 May;71(5):558-63.

Wang JC, To EW. Application of dermal substitute (Integra) to donor site defect of forehead flap. *Br J Plast Surg.* 2000 Jan;53(1):70-2.

King P. Artificial skin reduces nutritional requirements in a severely burned child. *Burns.* 2000;26:501-3.

Loss M, Wedler V, Kunzi W, Meuli-Simmen C, Meyer VE. Artificial skin, split thickness autograft and cultured autologous keratinocytes combined to treat a severe burn injury of 93% of TBSA. *Burns.* 2000;26:644-652.

Prystowsky JH, Nowygrod R, Marboe CC, Benvenisty AI, Ascherman JA, Todd GJ. Artificial skin (Integra Dermal Regeneration Template) for closure of lower extremity wounds. *Vascular Surgery.* 2000;34(6):557-567.

Fitton AR, Drew P, Dickson WA. The use of a bilaminar artificial skin substitute (Integra) in acute resurfacing of burns: an early experience. *Br J Plast Surg.* 2001 Apr;54(3):208-12.

Chan ES, Lam PK, Liew CT, Lau HC, Yen RS, King WW. A new technique to resurface wounds with composite biocompatible epidermal graft and artificial skin. *J Trauma.* 2001 Feb;50(2):358-62.

Chou T-D, Chen S-L, Lee T-W, Chen S-G, M.D.Cheng T-Y, Lee C-H, Chen T-M, Wang H-J. Reconstruction of Burn Scar of the Upper Extremities with Artificial Skin. *Plastic and Reconstructive Surgery.* 2001;108:378-384.

Moiesen NS, Staiano JJ, Ojeh NO, Thway Y, Frame JD. Reconstructive Surgery with a Dermal Regeneration Template: Clinical and Histologic Study. *Plastic and Reconstructive Surgery.* 2001;108:93-103.

Thomas WO, Rayburn S, Leblanc RT, Powell RW, Luterman A. Artificial Skin in the Treatment of a Large Congenital Nevus. *Southern Medical Journal.* 2001;94(3):325-328.

Prystowsky JH, Siegeal DM, Ascherman JA. Artificial skin for closure and healing of wounds created by skin cancer excision. *Dermatologic Surgery.* 2001;27:648-55.

Giovannini UM, Teot L. Aesthetic complex reconstruction of the lower leg: application of a dermal substitute (Integra) to an adipofascial flap. *Br J Plast Surg.* 2002 Mar;55(2):171-2.

Ryan CM, Schoenfeld DA, Malloy M, Schulz III JT, Sheridan RL, Tompkins RG. Use of Integra artificial skin is associated with decreased length of stay for severely injured adult burn survivors. *J Burn Care Rehabil.* 2002 Sep-Oct;23(5):311-7.

Fang P, Engrav LH, Gibran NS, Honari S, Kiriluk DB, Cole JK, Fleckman P, Heimbach DM, Bauer GJ, Matsumura H, Warner P. Dermatome Setting for Autografts to Cover Integra. *J Burn Care Rehabil.* 2002 Sep-Oct;23(5):327-32.

Análise de Pele Artificial

Tompkins RG, Hilton JF, Burke JF, Schoenfeld DA, Hegarty MT, Bondoc CC, Quinby WC Jr, Behringer GE, Ackroyd FW. Increased survival after massive thermal injuries in adults: preliminary report using artificial skin. *Crit Care Med.* 1989;17(8):734-40.

Jaksic T, Burke JF. The use of "artificial skin" for burns. *Annu Rev Med.* 1987;38(-HD-):107-17. Tompkins RG, Burke JF. Progress in burn treatment and the use of artificial skin. *World J Surg.* 1990 Nov-Dec;14(6):819-24.

Tompkins RG, Burke JF. Burn wound closure using permanent skin replacement materials. *World J Surg.* 1992 Jan-Feb;16(1):47-52.

Helvig EI. Dermal replacement: an update. *Semin Perioper Nurs.* 1997 Oct;6(4):233-5.

Cameron S. Changes in burn patient care. *Br J Theatre Nurs.* 1997 Aug;7(5):5-7.

Schultz JT, Tompkins RG, Burke JF. Artificial Skin. *Annu Rev Med.* 2000;51:231-244.

Orgill DP, Straus FH 2nd, Lee RC. The use of collagen-GAG membranes in reconstructive surgery. *Ann N Y Acad Sci.* 1999;888:233-48.

Winfrey ME, Cochran M, Hegarty MT. A new technology in burn therapy: Integra artificial skin. *Dimens Crit Care Nurs.* 1999;18(1):14-20.

Machens HG, Berger AC, Mailaender P. Bioartificial skin. *Cells Tissues Organs.* 2000;167(2-3):88-94.

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ANVISA: 10306840099
AWC-0003-V01

Referências

1. Singer AJ, Clark RAF. Cutaneous wound healing. *N Engl J Med.* 1999; 341:738–746.
2. Burke JF: Observations on the development and clinical use of artificial skin: An attempt to employ regeneration rather than scar formation in wound healing. *Jpn J Surg.* 1987;17(6) 431–438.
3. Heimbach D, Luterman A, Burke JF et al: Artificial dermis for major burns: A multi-center randomized clinical trial. *Ann Surg.* 1988 (Set);208(3):313–320.
4. Stern R, McPherson M, Longaker MT: Histologic study of artificial skin used in the treatment of full-thickness thermal injury. *J Burn Care Rehabil.* 1990;11:7–13.
5. Sheridan RL, Hegarty M, Tompkins RG, Burke JF: Artificial skin in massive burns—results to ten years. *Eur J Plastic Surgery.* 1994;17:91–93.
6. Burke JF, Yannas IV, Quinby WC, et al: Successful use of a physiologically acceptable artificial skin in the treatment of extensive burn injury. *Ann Surg.* 1981 (Out);194(4):413–428.
7. Berger A, Tanzella U, Machens HG, Liebau J, The use of Integra in primary burn wounds and unstable secondary scars. *Chirurg.* 2000; 71(5):558–563.
8. Tompkins RG, Hilton JF, Burke JF, Schoenfeld DA, Hegarty MT, Bondoc CC, Quinby WC, Behringer GE, Ackroyd FW. Increased survival after massive thermal injuries in adults: preliminary report using artificial skin. *Critical Care Medicine.* 1989; 17(8)734–740.
9. Data on file, Ethicon INC.
10. Besner G, Klamar J. Integra Artificial Skin as a Useful Adjunct in the Treatment of Purpura Fulminans. *J Burn Care Rehabil.* 1998;19:324–329.
11. Orgill DP, Butler CE, Regan JF. Behaviour of collagen-GAG matrices as dermal replacement in rodent and porcine models. *Wounds: A Comp of Clin. Research and Practice.* 1996;8(5):151–157.

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